

Innovative.
Simple.
Strong.



bag-aloneTM


part of the anchor **tissue retrieval system**



TRS070



TRS080



TRS090

anchor bag-alone™

features

- ✓ Offers safety & security
- ✓ Surgical defect can be kept small
- ✓ Virtual elimination of breaks & leaks
- ✓ Exceptional strength & quality
- ✓ Secure strong seal
- ✓ Rip-stop nylon

uses

	TRS 070	TRS 080	TRS 090
Cholecystectomy	✓	✓	
Appendectomy	✓	✓	
Oophorectomy	✓	✓	✓
Lobectomy	✓	✓	✓
Colectomy	✓	✓	✓
Prostatectomy	✓	✓	
Nephrectomy			✓
Splenectomy			✓
VATS			✓

dimensions

item	trocac size mm	mouth diameter cm/in	bag length cm/in	full capacity ml
TRS070	10	5.6/2.2	14.0/5.5	235
TRS080	10/12	7.6/3.0	14.0/5.5	300
TRS090	15	14.1/5.5	25.4/9.9	1550

ordering

item	description	UOM
TRS070	Tissue retrieval bag 070, Bag-Along™	10/BX
TRS080	Tissue retrieval bag 080, Bag-Along™	10/BX
TRS090	Tissue retrieval bag 090, Bag-Along™	5/BX

distributor



anchor™

about anchor

Since 1925, Anchor's mission has been to manufacture the finest surgical products. Expanding from the days of suture needles only, Anchor is now the only U.S. manufacturer to offer a full line of surgical specimen retrieval bags, both cannulated & non-cannulated. Anchor is dedicated to continuing innovation in specimen retrieval bag technology to meet the ever growing needs of surgeons worldwide.

ISO 13485: 2003 CMDCAS
ISO 13485: 2003
ISO 9001: 2008
CE 0297

MADE IN THE USA



Anchor Products Company, 52 Official Road, Addison, Illinois 60101 USA
Phone: 001-630-543-9124 | Fax: 001-630-543-9131 | info@anchorsurgical.com | www.anchorsurgical.com

The word 'Anchor', the symbol depicting the Caduceus entwined about an anchor comprising suture needles, and other distinctive terms, are trademarks of the Anchor Products Company, Addison, IL. The use of the color purple in conjunction with a surgical and medical apparatus or instrument is a registered trademark of Anchor Products Company, Addison, IL. Patents pending. All rights reserved.